



PressureGuard® Air Therapy Surface Matrix

Alternating Pressure/ Lateral Rotation		Microclimate Management ["Low Air Loss" (LAL)]		Bariatric	
	APM2	Protocol®	EasyAir® LR	Bariatric APM	Easy Air® XL
Therapy Modalities	Alternating Pressure	✓	✓	✓	✓
	Lateral Rotation	✓	✓	✓	✓
Power Outage Protection	Microclimate Management ("LAL")	✓	✓	✓	✓
	Powered Floatation	✓	✓	✓	✓
Control Unit Features	Can maintain inflation in power outage; (A) = Automatically	✓	✓ (A)	✓	✓ (A)
	Alarm: (V)isual, (A)djustable	V, A	V, A	V, A	V
Mattress Construction	Timed Auto Firm	✓	✓	✓	✓
	Geo-Matt® anti-shear top surface	✓	✓	✓	✓
	Safety Edge™ bolster design	✓	✓	✓	✓
	Heel Slope®, (F) = foam only heel section	✓	✓	✓	✓
	100% Latex free, 100% DHEP-free	✓	✓	✓	✓
Cover	Non-stretch coated nylon			✓	
	LifeSpan® bi-directional stretch with Shear Transfer Zone® design	✓			
	Zippered air flow coverlet with Air Diffused Matrix™ layer		✓		✓
Infection Control	Approved agents include bleach (1:10); phenolics; quaternary ammonia; hydrogen peroxide and alcohol.	✓	✓	✓	✓
General	Medicare Approved	E0277	E0277	E0277	E0277
	Therapeutic Weight Limit (lbs.)	350	400	750	1,000
	Meets Federal Fire Safety Regulations	✓	✓	✓	✓



All products shown manufactured by:
Span-America Medical Systems, Inc.